


- (en) HELMET FOR NON INVASIVE VENTILATION (N.I.V.)**
- (fr) CASQUE POUR VENTILATION NON INVASIVE (N.I.V.)**
- (de) HELM FÜR NICHT INVASIVE VENTILATIONSTHERAPIEN (N.I.V.)**
- (es) CASCO PARA VENTILACIÓN NO INVASIVA (N.I.V.)**
- (pt) CAPACETE PARA VENTILAÇÃO NÃO INVASIVA (V.N.I.)**
- (it) CASCO PER VENTILAZIONE NON INVASIVA (N.I.V.)**
- (nl) KAP VOOR NIET-INVASIEVE VENTILATIE (N.I.V.)**
- (no) PUSTEPOSE MEDUTVUNGENVENTILASJON (N.I.V. Non-Invasive Ventilation)**
- (fi) KYPÄRÄ NON-INVASIIVISTA VENTILAATIOTA (N.I.V.) VARTEN**
- (sv) HJÄLM FÖR NON-INVASIV VENTILATION (NIV)**
- (da) HJELM TIL IKKE-INVASIV VENTILATION**
- (el) KASKA GIA MH EPEMBATIKO MHCANIKO AERISMO (N.I.V.)**
- (pl) HE M DO WENTYLACJI NIEINWAZYJNEJ (N.I.V.)**
- (cs) MASKA PRO NEINVAZIVNÍ VENTILACI (N.I.V.)**
- (hu) SISAK NON-INVAZÍV (NIV) LÉLEGEZTETÉSHEZ**
- (sl) âELADA ZA NEINVAZIVNO ZRAâENJE (N.I.V.)**
- (et) KIIIVER MITTEINVASIIVSEKS VENTILATSIOONIKS (N.I.V.)**
- (ro) CASCĂ PENTRU VENTILAȚIE NEINVAZIVĂ (N.I.V.)**
- (sk) HELMA NA NEINVAZÍVNU VENTILÁCIU (N.I.V.)**
- (hr) KACIGA ZA NEINVAZIVNU VENTILACIJU (N.I.V.)**
- (tr) NON-INVAZIF VENTILASYON (N.I.V.) BAFILİ**
- (ar) خوذة للتنهوية غير الباضعة (N.I.V)**
- (he) קסדה להנשמה לא פולשנית (V.I.N)**
- (ua) Шолом для неінвазивної вентиляції (N.I.V.)**
- (sr) KACIGA ZA NEINVAZIVNU VENTILACIJU (N.I.V.)**

Instructions for use / Mode d'emploi / Gebrauchsanleitung / Instrucciones de uso / Instruções de utilização / Istruzioni per l'uso / Gebruiksaanwijzing / Brukerveiledning / Käyttöohjeet / Bruksanvisning / Brugsanvisning / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / Instrukcja użyci / Návod k použití / Használati útmutató / Navodilo za uporabo / Kasutusjuhend / Instrucțiuni de utilizare / Návod na použitie / Upute za uporabu / Sembollerin Anlamları / تعليمات للاستعمال / הוראות שימוש / Інструкція з використання / Uputstvo za Upotrebu



A1


A2

A2


A3


B1


B2


B3


B4


B6


C1


D1


D2


D3


D4



E2

E3

E4

(en)

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Intended use

Helmet for non-invasive ventilation in hospital environment.

Indications for use

Suitable for treatment of respiratory failure in hypoxicemic / hypercapnic patients. For other pathologies the treatment may only be performed under careful monitoring by expert operators.

The infant models are suitable for treatment of hypoxicemic patients.

Contraindications

- Coma
- uncooperative patient
- cardiac arrest
- hemodynamic instability
- recent oesophageal and gastro-surgical operations
- heavy bleeding of the upper digestive tract
- obstruction of the upper airways
- pneumothorax

Restrictions of use

1. The device must be used by qualified and trained medical/nursing staff; when used with a ventilator, the patient/ventilator interaction (trigger) may be complex and the device may therefore only be used by expert operators.
2. To be used for administration of air and oxygen.
3. When used for CPAP, use flow generators able to provide a total continuous air and oxygen flow of at least 40 L/min for adults and at least 30 L/min for paediatric patients and infants in order to ensure good lavage of the CO₂ exhaled. When used with a ventilator, make sure that it is able to deliver a flow during the inspiratory phase, sufficient to quickly remove the CO₂ from the inside the helmet.
4. If used properly, the device may be used continuously for a maximum of 7 days after which it must be replaced.
5. Monitor the clinical parameters of the patient. Warning: the models without anti-asphyxiation valve may only be used if monitoring of the clinical parameters is guaranteed with adequate systems equipped with alarms.
- 6.

PATIENTS	RECOMMENDED SIZE
ADULT - weighing > 30 kg	ADULT XS-S-M-L-XL-XXL, based on the neck circumference indicated on the label
PAEDIATRIC - To be used in paediatric age. - weighing approximately >12 kg (CPAP), >15 kg (NIV)	PAEDIATRIC, based on the neck circumference indicated on the label and the patient's age indicated on the product technical data sheet
INFANT CPAP - weighing approximately between 3 and 7 kg.	INFANT LOW, based on the approximate weight indicated on the product technical data sheet
INFANT CPAP - weighing approximately between 7 and 12 kg.	INFANT HIGH, based on the approximate weight indicated on the product technical data sheet

The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size.

Suitable choices will ensure a good seal around the patient's neck. The neck circumferences are approximate; the physician is advised to measure the neck circumference of the patient to be treated.

It is recommended to use the included measuring tape to choose the best size according to the patient's neck circumference.

Complications/side effects

Axillary pain (in patients with fragile skin and in prolonged therapies).

Claustrophobia.

CO₂ retention (rebreathing).

A. SETTING UP THE DEVICE

Remove the device from the package and model it with one hand in such a way that its cylindrical shape is restored. Before positioning the helmet on the patient open and close the access port (if present) in order to facilitate opening during use. Pre-check the valve. Pull and release the knob checking that the components slide smoothly.

1. Arrange the fastening systems securing them on the rear of the helmet.

2. When used with a flow distributor in CPAP

Connect the breathing circuit to the inlet connector. Connect the PEEP valve to the outlet connector. It is recommended to set a PEEP of at least

5cm/H₂O as prescribed by the physician. The integrated manometer can be useful to check if there is pressure in the helmet (the indication is however approximate). For the versions not equipped with an integrated manometer an external manometer can be connected to the pressure port or to an airtight access.

When used with a ventilator

Connect the breathing circuit to the inlet and the outlet connector.

3. Check that the connection is secure before proceeding with administration of the therapy. Activate the ventilation flow before continuing with the next positioning phase in order to reduce the time required to pressurise the helmet.

B. POSITIONING AND ACTIVATING THE DEVICE

1. Open out the collar using 4 hands, in such a way that the patient's head can pass through. To facilitate fitting the helmet, it is suggested that two people open the collar by gripping the flaps in diametrically opposite way and holding their thumbs on the rigid ring of the helmet.
2. Complete positioning of the fastening system on the front of the helmet.
3. Adjust the length of the straps so that the rigid ring is about 1 centimetre from the patient's shoulders.
4. Pressurise the system by pulling the knob of the anti-asphyxiation valve until the helmet reaches the desired internal pressure. The valve will go into patient safety and protection mode when, for any reason, the internal helmet pressure drops below 2 cm H₂O.
5. For proper functioning of the valve, do not obstruct the valve area and do not lock the action mechanism.
6. The helmet has airtight accesses to insert probes/catheters with diameters between 3.5 and 7 mm (ADULT and PAEDIATRIC sizes) and between 2.5 and 4 mm (INFANT size).

C. SYSTEMS TO IMPROVE COMFORT**1. Alternative fastening systems**

During extended therapy and with high operating pressure the patient may complain of axillary pain. Intersurgical can supply different accessories to relieve the patient:

- An abdominal/thoracic belt to be fixed as shown in Figure.
 - An elastic strap with loops to fix the helmet to the bed.
- To alleviate pain, hydrocolloid dressings may be applied between the patient's skin and the axillary protection.

2. Cushions / internal cuffs

- Some models have an internal inflatable cuff.
- The INFANT sizes have an anatomic cushion to rest the child's head on.

D. USING THE ACCESS PORT (when available)

1. Open the access port to the patient by turning the anti-asphyxiation valve anticlockwise.
2. Access the patient's face.
3. Close the access port to the patient by turning the anti-asphyxiation valve clockwise.
4. Pull the valve knob using two fingers until the helmet is pressurised.

E. HOW TO REMOVE THE HELMET

1. Remove any probes/catheters used during the therapy from the airtight accesses. Deflate the internal cuff (where present) by opening the clamp.
2. Detach the fastening systems.
3. Remove the helmet by opening out the collar with 4 hands.
4. Turn off the ventilation flow.

WARNINGS / PRECAUTIONS

Do not re-use. Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
2. Monitor the clinical parameters of the patient. Warning: models without anti-asphyxiation valve must only be used if monitoring of the clinical parameters is guaranteed with adequate systems equipped with alarms.
3. If the patient's conditions DETERIORATES or there is no improvement within the scheduled times, it is recommended to evaluate alternative ventilation Techniques.

- The efficacy of the therapy is considerably affected by the PEEP values set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient. Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.
- Caution: the device contains metal; do not use in departments where the presence of metal is a source of risk to the safety of the patient or third parties.
- During therapy, the patient may complain of annoying noise caused by the high flows utilised. This is due to the velocity and turbulence of the gas. The annoying noise can easily be eliminated or reduced by fitting one or two combined filters for breathing circuits on the helmet connectors - those normally used in hospital departments (not provided in this kit) - or having the patient wear ear plugs.
- When used with humidified ventilation gases, check that any condensate forming in the helmet does not cause the patient discomfort.
- The patient must REMOVE all accessories such as earrings, hair clips, combs and any other metal objects.
- The physician is responsible for deciding which ventilation therapy is most suited to the patient's pathology. The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. The neck circumferences (expressed in centimetres) and the helmet size is indicated on each device. Suitable choices will ENSURE A GOOD SEAL around the patient's neck.
- Use ONLY ON UNBROKEN skin.
- The device can be used for maximum 7 days. Depending on the clinical picture of the patient, the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
- Non-sterile. Do not sterilize.
- Do not use on another patient.
- Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°/+50°C).
- To facilitate modelling/shaping the helmet, it is recommended to keep it at room temperature (about 20-25°C) for a few hours before use.
- Dispose of the materials immediately after use CONFORMING TO the current laws and regulations.
- The cuff inflation pump is not included in the kit. Do not use pressurised gas sources to inflate the integrated cuff.
- Activating the anti-asphyxiation valve in the event of ventilation interruption allows air exchange with the outside limiting CO₂ rebreathing. The anti-asphyxiation valve does not substitute the ventilation support in the event of accidental ventilation interruption. **The anti-asphyxiation valve does not exempt the nursing staff from implementing appropriate monitoring and supervision.**
- Caution: Use of the helmet does not exempt operators from the obligation to use adequate Personal Protection Devices in accordance with the hospital procedures.
- Non invasive ventilation shall be interrupted if there is one of the following conditions:
 - worsening state of consciousness and respiratory distress
 - loss of airway protection
 - unchanged PaCO₂ (during two subsequent EGA performed after max 1 hour)
 - persistency of severe hypoxemia
 - severe and uncontrollable hemodynamic instability
 - patient / ventilator synchronization problems
 - uncontrollable secretions
 - interface intolerance

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy