VENTUKIT

- it** CASCO PER N.I.V. (CPAP) CON SISTEMA DI MISCELAZIONE ARIA/OSSIGENO
- en** HELMET FOR N.I.V. (CPAP) WITH AIR/OXYGEN BLENDING SYSTEM
- fr** CASQUE POUR N.I.V. (CPAP) AVEC SYSTEME DE MELANGE AIR/OXYGENE
- de** HELM FÜR N.I.V. (CPAP) MIT SYSTEM ZUR LUFT-SAUERSTOFF-MISCHUNG
- es** CASCO PARA N.I.V. (CPAP) CON SISTEMA DE MEZCLADO DE AIRE/OXÍGENO
- pt** CAPACETE PARA V,N.I. (CPAP) COM SISTEMA DE MISTURA AR/OXIGÉNIO
- nl** KAP VOOR N.I.V. (CPAP) MET LUCHT/ ZUURSTOFMENGSTEEEM
- sv** HJÄLM FÖR NIV (CPAP) MED BLANDNINGSSYSTEM LUFT/SYRGAS
- el** ΚΑΣΚΑ ΜΕ ΣΥΣΤΗΜΑ ΑΝΑΜΙΞΗΣ ΑΕΡΑ/ΟΞΥΓΟΝΟΥ
- sl** ČELADA S SISTEMOM ZA MEŠANO DOVAJANJE ZRAKA/KISIKA
- sr** KAČIGA ZA N.I.V. (CPAP) SA SISTEMOM ZA MEŠANJE VAZDUHA / KISEONIKA
- sk** HELMA PRE N.I.V. (CPAP) SO SYSTÉMOM NA MIEŠANIE VZDUCHU/KYSLÍKA
- tr** NIV (CPAP) için hava/oksijen karışım sistemli başlık
- ua** ШОЛОМ ДЛЯ НІВ (CPAP) зі вбудованою системою змішування повітря з киснем.

Istruzioni per l'uso / Instructions for use / Mode d'emploi / Gebrauchsanleitung / Instrucciones de uso / Instruções de utilização / Gebruiksaanwijzing / Bruksanvisning / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / Navodilo za uporabo / UPUTSTVO ZA UPOTREBU / Návod na pouzitie / Kullanma Talimatları / ІНСТРУКЦІЇ З ВИКОРИСТАННЯ

LEGENDA SIMBOLI - SYMBOL LEGEND - LEGENDE DES SYMBOLES - ZEICHENERKLÄRUNG - LEYENDA DE LOS SIMBOLOS - LEGENDA DOS SÍMBOLOS - VERKLARING DER TEKENS - TECKENFÖRKLARING - ΥΠΟΜΗΜΑ ΣΥΜΒΟΛΩΝ - LEGENDA SIMBOLI - LEGENDA SIMBOLA - LEGENDA SYMBOLOV - Sembollerin Anlamları - УМОВНИ ПОЗНАЧЕННЯ

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KIT CONTENTS

- Non-invasive ventilation hood (CPAP) with integrated air/oxygen mixing system.
- Oxygen feed tubes
- PEEP valve and pressure gauge
- Anti-noise earplugs
- Instructions for use.

Note: the flowmeters are not included in the kit.

Intended use

Ventukit is a device intended for non-invasive CPAP ventilation (Continuous Positive Airway Pressure) in a hospital environment when an FiO_2 variable between 40% and 100% needs to be supplied.

Indications for use

Suitable for treatment of patients with acute hypoxemic respiratory failure. For other pathologies the treatment may only be performed under careful monitoring by expert operators.

Contraindications

- Coma
- Uncooperative patient
- Cardiac arrest
- Haemodynamic instability
- Recent oesophageal and gastro-surgical operations
- Heavy bleeding of the upper digestive tract
- Obstruction of the upper airways
- Pneumothorax

Restrictions of use

1. The device must be used by qualified and trained medical/nursing staff.
2. To be used for administration of air and oxygen.
3. For pure oxygen treatment, connect both tubes to oxygen sources and feed the device with O_2 flows of not less than 40 l/min deriving from the sum of the two lines.
4. If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.
5. Monitor the clinical parameters of the patient. Warning: models without an anti-asphyxiation valve may be used only if monitoring of the clinical parameters can be guaranteed with adequate systems equipped with alarm.
6. To be used solely with oxygen flowmeters conform to the standards and capable of delivering the flows indicated in the table. If fed with an O_2 source (only line A), Ventukit can reach an FiO_2 of 40%. To achieve a higher FiO_2 use two O_2 sources at the same time (line A and line B). The use of flow meters with an output of less than the maximum required in the table on line B, will reduce the performance of the device.
7. Use flowmeters with a 7mm diameter hose tail.
- 8.

PATIENTS	RECOMMENDED SIZE
ADULTS - weighing > 30 kg	ADULT XS-S-M-L-XL-XXL, based on the neck circumference indicated on the label

The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. An appropriate choice will ensure good hermetic seal around the patient's neck. The neck circumferences are approximate; the physician is advised to measure the neck circumference of the patient to be treated. **It is recommended to use the included measuring tape to choose the best size according to the patient's neck circumference.**

Where and by whom can Ventukit be used

Ventukit may only be used by qualified medical/nursing staff in a hospital environment. The therapy is started simply by connecting Ventukit to the oxygen source available in hospital wards (wall-mounted distribution systems or cylinders), for example: emergency medicine, pneumology, haematology, oncology, post-surgery intensive care, ambulance, emergency room, etc.

Complications and side effects

Episodes of armpit pain may occur in patients that have a sensitive skin and during extended therapies at high PEEP pressures, which can be avoided by applying thin strips of hydrocolloid or by alternating the fastening systems such as chest straps or elastic bands for fastening to the bed. Some patients may

suffer episodes of claustrophobia. The helmet is equipped with a quick-access opening so that the operator can easily access the patient's face. Improper use of the device resulting in flows of less than 40 l/min may cause an increase of CO_2 inside the helmet. Ventukit, if used in compliance with these instructions, guarantees a total flow to the patient of more than 40 l/min thus obtaining proper ventilation gas exchange. During extended CPAP therapies, the noise inside the helmet might create some discomfort to the patient. To prevent this, Ventukit is equipped with an integrated silencer that can reduce the noise generated by the Venturimeter by approximately 10 dB. The kit also contains earplugs which provide considerable acoustic comfort for the patient.


Fig. 1


Fig. 2

How to use the VENTUKIT

The therapy is started simply by connecting VENTUKIT to the available oxygen source.

The Ventukit tubes A and B then need to be properly connected to the respective flowmeters (Ref. Figs. 1 and 2) connected to the oxygen source in the ward appropriately pressurised following the instructions provided by the manufacturer of the flowmeters. Ventukit is provided with a performance table of the system's operating conditions.

Example: to administer 50% FiO_2 + a PEEP of 7.5 cm H_2O set the flowmeter A to 8 l/min and the flowmeter B to 9 l/min, thus obtaining a flow of about 47 l/min to the patient.

The performances of the Venturimeters are by their nature affected by the resistance variations induced by the PEEP valves. It follows that the flow and concentration values may undergo some variations. For a more correct flow control at different PEEP values, it is recommended to refer to the present document. Please note that flowmeter A is used for functioning of the Venturimeter: **+ flow in A = + flow to the patient**

Flowmeter B, on the other hand, allows increasing the concentration of oxygen: **+ flow in B = + FiO_2**

Precautions for use of the flowmeters

After using the device, it is essential to turn the flowmeter stopcocks until they are fully closed so that gas delivery is interrupted only after having removed Ventukit from the patient.

For use of the flowmeters, always refer to the individual instructions for use contained in the product package.

100%	To deliver 100% FiO_2 , occlude the air inlet of the Venturi with the blue collar. Connect both A and B lines to an oxygen source and provide a total flow of at least 40 L/min.
-------------	--

Refer to the table attached to the device

Warnings

It is recommended to use the settings in the table provided by Intersurgical. These settings assure sufficient flows for proper CPAP therapy with helmet. Higher flow settings would generate unnecessary oxygen consumption and excessive system noise. Should parameters different from those listed in the table be set, it is recommended to check the therapy conditions with the appropriate instruments (oximeters, flowmeters, etc.).


HOW TO ADMINISTER 100% OXYGEN TO THE PATIENT

Fig. 3

If 100% FiO_2 needs to be supplied, it is essential to deactivate operation of the Venturimeter.

This is achieved by sliding the blue ring nut "G" onto the air inlet of the Venturimeter (ref. Figs. 3 and 4)

The total flow to the patient will be the sum of the flows indicated by

the flowmeters A and B.

Fig. 4

To restore operation of the Venturimeter and thus reduce the FiO_2 delivered, slide back the ring nut as shown in Fig. 3 restoring the therapy conditions according to the present document.

Warnings

Always check that the ventilation circuit is perfectly connected. In the event that the tubes disconnect, patient ventilation will fail. It is recommended to frequently monitor the state of the patient during the therapy. To assure that the FiO_2 delivered to the patient is accurately measured, connect an oximeter to the monitoring port and run random tests.

A. SETTING UP THE DEVICE

Remove the device from the package and model it with one hand in such a way that its cylindrical shape is restored. Before putting the helmet on the patient, open and close the access port (if present) in order to facilitate opening during use. Pre-check the valve. Pull and release the knob checking that the components slide smoothly.

1. Arrange the fastening systems securing them on the rear of the helmet.
2. Connect the lines A and B to flow as explained above. It is recommended to set a PEEP of at least 5cm/ H_2O as prescribed by the physician. The integrated manometer can be useful to check if there is pressure in the helmet (the indication is however approximate). For the versions not equipped with an integrated manometer an external one can be connected to the pressure port or to an airtight access.
3. Check that the airtight accesses are correctly positioned. Set up the parameters as for the included table, according to physician's prescription. Activate the ventilation flow before continuing with the next positioning phase in order to reduce the time required to pressurise the helmet.
4. During extended CPAP therapies, the noise inside the helmet might create some discomfort to the patient. To prevent this, Ventukit is equipped with an integrated silencer that can reduce the noise generated by the Venturimeter by approximately 10 dB. The kit also contains earplugs which provide considerable acoustic comfort for the patient.

B. POSITIONING AND ACTIVATING THE DEVICE

1. Open out the collar in such a way that the patient's head can pass through. To facilitate fitting the helmet it is suggested that two operators open out the collar by gripping the flaps in diametrically opposite way and holding their thumbs on the rigid ring of the helmet.
2. Complete positioning of the fastening system on the front of the helmet.
3. Adjust the length of the straps so that the rigid ring is about 1 centimetre from the patient's shoulders.

For models with anti-asphyxiation valve

4. Pressurise the system by pulling the knob of the anti-asphyxiation valve until the helmet reaches the desired internal pressure. The valve will revert to a safe state to protect the patient when, for any reason, the internal helmet pressure drops to below 2 cm H_2O .
5. For proper functioning of the valve, do not obstruct the valve area and do not lock the action mechanism.
6. The helmet has airtight accesses to insert probes/catheters with diameters between 3.5 and 7 mm.

C. SYSTEMS TO IMPROVE COMFORT

Alternative fastening systems

1. During extended therapy and with high operating pressures the patient may complain of armpit pain. Intersurgical can supply different accessories to give the patient relief:
 - An abdominal/chest belt to be fixed as shown in the figure.
 - An elastic strap with slots to fix the helmet to the bed.
 To alleviate pain, hydrocolloid dressings may be applied between the patient's skin and the axillary protection.

D. USING THE ACCESS PORT (if available)

1. Open the access port to the patient by turning the anti-asphyxiation valve anticlockwise.
2. Access the patient's face.
3. Close the access port to the patient by turning the anti-asphyxiation valve clockwise.
4. Pull the valve knob using two fingers until the helmet is pressurised.

E. HOW TO REMOVE THE HELMET

1. Remove any probes/ catheters used during the therapy from the airtight accesses.
2. Detach the fastening systems.
3. Remove the helmet by opening out the collar (two operators).
4. Turn off the ventilation flow.

WARNINGS / PRECAUTIONS

⚠ Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

⚠ Caution: this product contains phthalate. There is evidence suggesting that phthalate exposure from medical procedures may harm children, pregnant women or breastfeeding women. Use of this product does not represent a high exposure risk to phthalate, however, as a precautionary measure its use should be limited to essential procedures.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
2. Monitor the clinical parameters of the patient. Warning: models without an anti-asphyxiation valve may be used only if monitoring of the clinical parameters can be guaranteed with adequate systems equipped with alarm.
3. If noticing that the patient's conditions deteriorate or there is no improvement within the scheduled times, it is recommended to evaluate alternative ventilation techniques.
4. The efficacy of the therapy is considerably affected by the PEEP value set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient. Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.
5. Caution: The device contains metal; do not use in departments where the presence of metal is a source of risk for the safety of the patient or third parties.
6. During therapy, the patient may complain of annoying noise caused by the high flows utilised. This is due to the velocity and turbulence of the gas. The annoying noise can easily be reduced by fitting one or two combined filters for breathing circuits on the hood connectors (normally found in hospital departments and not provided in this kit) or having the patient wear sound-absorbing ear plugs (provided in this kit).
7. When used with humidified ventilation gas, check that any condensate that forms does not discomfort the patient.
8. The patient must remove all accessories such as earrings, hair clips, combs and any other metal objects.
9. The physician is responsible for deciding which ventilation therapy is most suited to the patient's pathology. The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. The neck circumference (expressed in centimetres) and the helmet size is indicated on each device. An appropriate choice will ensure good hermetic seal around the patient's neck.
10. Use only on unbroken skin.
11. The device can be used for maximum 7 days. Depending on the clinical picture of the patient the medical/nursing staff are responsible for defining the need for more frequent replacement of the device.
12. Non-sterile. Do not sterilize.
13. Do not use on another patient.
14. Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).
15. To facilitate modelling/shaping the helmet, it is recommended to keep it at room temperature (about 20-25°C) for a few hours before use.
16. Dispose of the materials immediately after use in compliance with the legislation in force.
17. Activating the anti-asphyxiation valve in the event of ventilation interruption allows air exchange with the outside limiting CO_2 rebreathing. The anti-asphyxiation valve does not substitute the ventilation support in the event of accidental ventilation interruption. **The anti-asphyxiation valve does not exempt the nursing staff from implementing appropriate monitoring and supervision.**
18. Caution: Use of the helmet does not exempt operators from the obligation to use adequate Personal Protection Devices in accordance with the hospital procedures.
19. Non invasive ventilation shall be interrupted if there is one of the following conditions:
 - worsening state of consciousness and respiratory distress
 - loss of airway protection
 - unchanged PaCO_2 (during two subsequent EGA performed after max 1 hour)
 - persistence of severe hypoxemia
 - severe and uncontrollable hemodynamic instability
 - patient / ventilator synchronization problems
 - uncontrollable secretions
 - interface intolerance

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy