
CASTAR NEXT

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Intended use

Helmet for non-invasive ventilation in a hospital environment.

Indications for use

Suitable for treatment of respiratory failure in hypoxemic and/or hypercapnia patients. For other pathologies, the treatment may only be performed after careful clinical evaluation by expert operators.

Contraindications

- Coma
- Uncooperative patient
- Cardiac arrest
- Hemodynamic instability
- Recent oesophageal and gastro-surgical operations
- Heavy bleeding of the upper digestive tract
- Obstruction of the upper airways
- Pneumothorax
- It is inadvisable to use the device on patients with critical carotid stenosis

Restrictions of use

1. The device must be used by qualified and trained medical/nursing staff. The patient/ventilator interaction (trigger) may be complex and the device may therefore only be used by expert operators.
2. To be used for administration of air and oxygen.
3. When used for CPAP, use flow generators able to provide a total continuous air and oxygen flow of at least 40 L/min in order to ensure good lavage of the CO₂ exhaled. When used with a ventilator, make sure that it is able to deliver a flow, during the inspiratory phase, sufficient to quickly remove the CO₂ from inside the helmet.
4. If properly used, the device may be used continuously for 7 days after which it must be replaced.
5. Monitor the clinical parameters of the patient.
- 6.

PATIENTS	RECOMMENDED SIZE
ADULT - approximate weight > 30 kg	ADULT XS-S-M-L-XL-XXL, based on the neck circumference indicated on the label

The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. A suitable choice will ensure good seal around the patient's neck.

The neck circumferences are approximate; the physician is advised to measure the neck circumference of the patient to be treated. It is recommended to use the

included measuring tape to choose the best size according to the patient's neck circumference.

Complications/side effects

- Claustrophobia.
- CO₂ retention (rebreathing).

A. SETTING UP THE DEVICE

Remove the device from the package and model it with one hand in such a way that its cylindrical shape is restored.

Open the access port before placing the helmet on the patient. Pre-check the valve. Pull and release the knob checking that the components slide smoothly. Check the cushion by inflating it with the pump supplied. Deflate the cushion before positioning the helmet on the patient.

1. When used with a flow driver in CPAP

Connect the breathing circuit to the inlet connector **(A)**. Connect the PEEP valve to the outlet connector **(B)**. It is recommended to set a PEEP of at least 5cm/H₂O as prescribed by the physician. The integrated manometer can be useful to check if there is pressure in the helmet (the indication is however approximate). For the versions not equipped with an integrated manometer an external one can be connected to the pressure port or to an airtight access.

When used with a ventilator

Connect the breathing circuit to the inlet connector **(A)** and to the outlet connector **(B)**.

2. Check that the airtight accesses are correctly positioned. Check that the connection is secure before proceeding with administration of the ventilation gases. Activate the ventilation flow before continuing with the next positioning phase in order to reduce the time required to pressurise the helmet.

B. POSITIONING AND ACTIVATING THE DEVICE

1. Open out the collar **(II)** using 4 hands in such a way that the patient's head can pass through. To facilitate fitting the helmet, it is suggested that two people open the collar by gripping the flaps in diametrically opposite way and holding their thumbs on the rigid ring of the helmet.
2. The helmet has airtight accesses to insert probes/catheters with diameters between 3.5 and 7 mm.
3. **Access next to the port:** Remove the access from its rigid seat and open it. Position the probe/catheter and refit the access in its rigid seat. Alternatively, remove only the cap and use the central hole. Use this access for NGT if in use.

4. Position the inflatable cushion **(H)** around the patient's neck and hook the fasteners **(F2)** onto the knobs **(E)**.
- 5a. Take the rigid annular body **(G)**, open it and position it around the patient's neck following the indications on the body.
- 5b. Close the rigid annular body on the back of the patient's neck.
6. Using the fastening tabs **(F)** hook the rigid annular body **(G)** onto the knobs **(E)** positioned on the helmet ring.
7. Check that the cushion inflation line **(M)** runs underneath the rigid annular body **(G)** and fit it into place in the tube seat **(L)** in the lower part of the rigid annular body.
8. Use a pump **(N)** to inflate the cushion **(H)** through the inflation line **(M)** and close the clamp. The cushion serves to improve patient comfort and to hold the helmet in place during the treatment under pressure. Make sure that the elastic collar **(I)** juts out of the central hole of the cushion **(H)**.
9. Close the access port **(C)**. Pressurise the system by pulling the knob of the anti-asphyxiation valve until the helmet reaches the desired internal pressure. The valve reverts to a safe state to protect the patient when, for any reason, the internal helmet pressure drops to below approximately 2 cmH₂O.
10. For proper functioning of the valve, do not obstruct the valve area and do not lock the action mechanism.

C. USING THE ACCESS PORT (when available)

1. Open the access port to the patient by twisting off the anti-asphyxiation valve anticlockwise.
2. Access the patient's face.
3. Close the access port to the patient by twisting on the anti-asphyxiation valve clockwise.
4. Pull the valve knob using two fingers until the helmet is pressurised.

D. HOW TO REMOVE THE HELMET

1. Open the access port.
2. Remove any probes/catheters used during the therapy from the airtight accesses **(D)**.
3. Remove the cushion inflation line **(M)** from its seat **(L)**. Undo the fastening tabs **(F)** of the rigid annular body and the cushion fasteners **(F2)** from the front metal connectors **(E)** on the helmet.
4. Remove the rigid annular body **(G)** and move the cushion **(H)** to the back of the patient's neck. Remove the helmet by opening out the collar with 4 hands.
5. Turn off the ventilation flow.

WARNINGS / PRECAUTIONS

 Do not re-use. Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
2. Monitor the clinical parameters of the patient. Warning: models without anti-asphyxiation valve must only be used if monitoring of the clinical parameters is guaranteed with adequate systems equipped with alarms.
3. If the patient's conditions DETERIORATES or there is no improvement within the scheduled times, it is recommended to evaluate alternative ventilation Techniques.
4. The efficacy of the therapy is considerably affected by the PEEP values set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient. Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.
5. Caution: the device contains metal; do not use in departments where the presence of metal is a source of risk to the safety of the patient or third parties.
6. During therapy, the patient may complain of annoying noise caused by the high flows utilised. This is due to the velocity and turbulence of the gas. The annoying noise can easily be eliminated or reduced by fitting one or two combined filters for breathing circuits on the hood connectors - those normally used in hospital departments (not provided in this kit) - or having the patient wear ear plugs.
7. When used with humidified ventilation gases, check that any condensate forming in the hood does not cause the patient discomfort.
8. The patient must REMOVE all accessories such as, earrings, hair clips, combs and any other metal objects.
9. The physician is responsible for deciding which ventilation therapy is most suited to the patient's pathology. The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. The neck circumferences (expressed in centimetres) and the hood size is indicated on each device. Suitable choices

will ENSURE A GOOD SEAL around the patient's neck.

10. Use ONLY ON UNBROKEN skin.
11. The device can be used for maximum 7 days. Depending on the clinical picture of the patient, the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
12. Non-sterile. Do not sterilize.
13. Do not use on another patient.
14. Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°/+50°C).
15. To facilitate modeling/shaping the hood, it is recommended to keep it at room temperature (about 20-25°C) for a few hours before use.
16. Dispose of the materials immediately after use CONFORMING TO the current laws and regulations.
17. Do not use pressurised gas sources to inflate the integrated cuff.
18. Activating the anti-asphyxiation valve in the event of ventilation interruption allows air exchange with the outside limiting CO₂ rebreathing. The anti-asphyxiation valve does not substitute the ventilation support in the event of accidental ventilation interruption. **The anti-asphyxiation valve does not exempt the nursing staff from implementing appropriate monitoring and supervision.**
19. Caution: Use of the helmet does not exempt operators from the obligation to use adequate Personal Protection Devices in accordance with the hospital procedures.
20. Non invasive ventilation shall be interrupted if there is one of the following conditions:
 - worsening state of consciousness and respiratory distress
 - loss of airway protection
 - unchanged PaCO₂ (during two subsequent EGA performed after max 1 hour)
 - persistence of severe hypoxemia
 - severe and uncontrollable hemodynamic instability
 - patient / ventilator synchronization problems
 - uncontrollable secretions
 - interface intolerance.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy