VENTUKIT NEXT

- it** CASCO PER N.I.V. (CPAP) CON SISTEMA DI MISCELAZIONE ARIA/OSSIGENO
- en** HELMET FOR N.I.V. (CPAP) WITH AIR/OXYGEN BLENDING SYSTEM
- fr** CASQUE POUR N.I.V. (CPAP) AVEC SYSTEME DE MELANGE AIR/OXYGENE
- de** HELM FÜR N.I.V. (CPAP) MIT SYSTEM ZUR LUFT-SAUERSTOFF-MISCHUNG
- es** CASCO PARA N.I.V. (CPAP) CON SISTEMA DE MEZCLADO DE AIRE/OXÍGENO
- nl** KAP VOOR N.I.V. (CPAP) MET LUCHT/ZUURSTOFMENGSTEEEM

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Instructions for use

Mode d'emploi

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
LEGENDA SIMBOLI - SYMBOL LEGEND - LEGENDE DES SYMBOLES ZEICHENERKLÄRUNG - LEYENDA DE LOS SIMBOLOS - VERKLARING DER TEKENS

MOD	Modello/taglia - Model/size - Modèle/Taille - Modell/Größe - Modelo/Tamaño - Model/maat		
REF	Codice - Code - Code - Code - Código - Code		
LOT	Lotto - Batch - Lot - Charge - Lote - Charge		Fragile - Handle with care - Fragile - Vorsicht zerbrechlich - Frágil - Voorzichtig behandelen
	Scadenza - Expiry - Date de péremption - Verfalldatum - Caducidad - Vervaldatum		Tenere al riparo dai raggi solari - Keep away from sunlight - Conserver à l'abri de la lumière - Vor Sonnenlicht schützen - Proteger de la luz del sol - Uit het zonlicht houden
QTY	Quantità - Quantity - Quantité - Menge - Cantidad - Aantal		Tenere al riparo dalla pioggia - Keep away from rain - Conserver à l'abri de l'humidité - Vor Regen schützen - Proteger de la lluvia - Uit de regen houden
	Non riutilizzare - Do not re-use - Ne pas utiliser - Nicht verwenden - No utilizar - Niet gebruiken		Privo di lattice - latex free - Ne contient pas de latex - Latexfrei - Sin látex - Bevat geen latex
	Attenzione - Caution - Mise en garde - Vorsicht! - Precaución - Waarschuwing		<u>Se presente su etichetta/If present on the label/si elle est présente sur l'étiquette/falls vorhanden auf dem Etikett/si está presente en la etiqueta/indien aanwezig op het etiket</u> Avvertenza: questo prodotto contiene ftalato - Caution: this product contains phthalate - Attention: ce produit contient des phtalates - Achtung: Dieses Produkt enthält Phthalate - Precaución: El producto contiene ftalatos - Waarschuwing: Er is ftalaat in dit product.
	Vedere le istruzioni d'uso allegate al dispositivo - See instructions for use attached to the device - Voir le mode d'emploi joint au dispositif - Siehe der Vorrichtung beiliegende Gebrauchsanleitung - Véase las instrucciones de uso que se adjuntan con el dispositivo - Zie de bij het hulpmiddel geleverde gebruiksaanwijzingen		
	Fabbricante - Manufacturer - Fabricant - Hersteller - Fabricante - Fabrikant		<u>Se presente su etichetta/If present on the label/si elle est présente sur l'étiquette/falls vorhanden auf dem Etikett/si está presente en la etiqueta/indien aanwezig op het etiket</u> Privo di ftalati - phtalates free - Ne contient pas de phtalates - Phthalate frei - Sin ftalatos - Bevat geen ftalaten
	Non usare se la confezione è danneggiata - Do not use if the package is damaged - Ne pas utiliser si le conditionnement est endommagé - Nicht verwenden, falls die Verpackung beschädigt ist - No usar si el envase está dañado - Niet gebruiken als de verpakking aangebroken is		
	Non sterile - Non-sterile - Non stérile - Nicht steril - No estéril - Niet steriel	MD	Dispositivo medico - Medical Device - Dispositif médical - Medizinprodukt - Dispositivo médico - Medisch hulpmiddel
	Limitazioni della temperatura - Temperature limitations - Seuls de temperature - Temperatureinschränkungen - Límites de temperatura - Temperatuurbependingen	UDI	Identificatore dispositivo univoco - Unique Device Identifier - Identifiant unique de l'appareil - Einmalige Produktkennung - Identificador único de dispositivo - Eenduidige ID van het hulpmiddel
	Circonferenza collo - Neck circumference - Circonférence du cou - Halsumfang - Circunferencia del cuello - Halsomtrek		
	Non aprire l'imballaggio con un coltello - No not open packaging using a knife - Ne pas ouvrir le conditionnement avec un couteau - Verpackung nicht mit Messer öffnen - No abrir el embalaje con un cuchillo - De verpakking niet openen met een mes		
	La marcatura CE include il numero identificativo TUV Rheinland Italia (Organismo Notificato). Prodotto conforme ai requisiti MDR (UE) 2017/745 - The CE marking includes the TUV Rheinland Italia (notified body) identification number. The product conforms to the requirements set out in the EU MDR 2017/745 - Le marquage CE inclut le numéro d'identification TUV Rheinland Italia (Organisme notifié). Produit conforme aux exigences du MDR EU 2017/745 - Zum CE-Zeichen gehört die TUV Rheinland Italia-Identifikationsnummer (benannte Stelle). Das Produkt entspricht den Anforderungen der Verordnung (EU) 2017/745 über Medizinprodukte - El marcado CE incluye el número de identificación TUV Rheinland Italia (Órgano Notificado). El producto		


A1

A2


A3


B1


B2


B3


B4


B5a


B5b


B6


B7


B8


B9


B9


C1


C2


C3


C4


D1


D2


D3

D4

D5

it

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en

- A** Inspiration connector with FiO₂ port
- B** PEEP valve connector
- C** Anti-asphyxiation valve / Inspection port (if present)
- D** Airtight accesses for probes or catheters
- E** Manometer
- F** 100% FiO₂ control ring nut
- G** Venturimeter
- H** Flowmeter connection tubes
- I** Knobs for fastening the cushion and the rigid annular body (also on the rear of the helmet)
- L** Rigid annular body fastening tabs
- M** Inflatable cushion fasteners
- N** Rigid annular body
- O** Inflatable cushion
- P** Collar
- Q** Cushion inflation tube seat
- R** Cushion inflation line
- S** Pump (provided)

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KIT CONTENTS

- Non-invasive ventilation hood (CPAP) with integrated air/oxygen mixing system.
- Oxygen feed tubes
- PEEP valve and pressure gauge
- Anti-noise earplugs
- Pump
- Instructions for use.

Note: the flowmeters are not included in the kit.

Intended use

Ventukit NEXT is a device intended for non-invasive CPAP ventilation (Continuous Positive Airway Pressure) in a hospital environment when an FI_{O_2} variable between 40% and 100% needs to be supplied.

Indications for use

Suitable for treatment of respiratory failure in hypoxemic patients. For other pathologies the treatment may only be performed under careful monitoring by expert operators.

Contraindications

- Coma
- Uncooperative patient
- Cardiac arrest
- Haemodynamic instability
- Recent oesophageal and gastro-surgical operations
- Heavy bleeding of the upper digestive tract
- Obstruction of the upper airways
- Pneumothorax
- It is unadvisable to use the device on patients with critical carotid stenosis.

Restrictions of use

1. The device must be used by qualified and trained medical/nursing staff.
2. To be used for administration of air and oxygen.
3. For pure oxygen treatment, connect both tubes to oxygen sources and feed the device with O_2 flows of not less than 40 l/min deriving from the sum of the two lines.
4. If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.
5. Monitor the clinical parameters of the patient. Warning: models without an anti-asphyxiation valve may be used only if monitoring of the clinical parameters can be guaranteed with adequate systems equipped with alarm.
6. To be used solely with oxygen flowmeters conform to the standards and capable of delivering the flows indicated in the table. If fed with an O_2 source (only line A), Ventukit Next can reach an FI_{O_2} of 40%. To achieve a higher FI_{O_2} use two O_2 sources at the same time (line A and line B). The use of flow meters with an output of less than the maximum required in the table on line B, will reduce the performance of the device.
7. Use flowmeters with a 7mm diameter hose tail.
- 8.

PATIENTS	RECOMMENDED SIZE
ADULT - weighing > 30 Kg	ADULT XS-S-M-L-XL-XXL, based on the neck circumference indicated on the label

The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. Suitable choices will ensure a good seal around the patient's neck.

The neck circumferences are approximate; the physician is advised to measure the neck circumference of the patient to be treated.

It is recommended to use the included measuring tape to choose the best size according to the patient's neck circumference.

Where and by whom can Ventukit NEXT be used

Ventukit NEXT may only be used by qualified medical/nursing staff in a hospital environment. The therapy is started simply by connecting Ventukit NEXT to the oxygen source available in hospital wards (wall-mounted distribution systems or cylinders), for example: emergency medicine, pneumology, haematology, oncology, post-surgery intensive care, ambulance, emergency room, etc.

Complications and side effects

Some patients may suffer episodes of claustrophobia. The helmet is equipped with a quick-access opening so that the operator can easily access the patient's face. Improper use of the device resulting in flows of less than 40 l/min may cause an increase of CO_2 inside the helmet. Ventukit NEXT, if used in compliance with these instructions, guarantees a total flow to the patient of more than 40 l/min thus obtaining proper ventilation gas exchange. During extended CPAP therapies, the noise inside the helmet might create some discomfort to the patient. To prevent this, Ventukit NEXT is equipped with an integrated silencer that can reduce the noise generated by the Venturimeter by approximately 10 dB. The kit also contains

earplugs which provide considerable acoustic comfort for the patient.

How to use the VENTUKIT NEXT

Fig. 1


Fig. 2

The therapy is started simply by connecting VENTUKIT to the available oxygen source.

The Ventukit NEXT tubes A and B then need to be properly connected to the respective flowmeters (Ref. Figs. 1 and 2) connected to the oxygen source in the ward appropriately pressurized following the instructions provided by the manufacturer of the flowmeters. Ventukit NEXT is provided with a performance table of the system's operating conditions. Example: to administer 50% FI_{O_2} + a PEEP of 7.5 cm H_2O set the flowmeter A to 8 l/min and the flowmeter B to 9 l/min, thus obtaining a flow of about 47 l/min to the patient.

The performances of the Venturimeters are by their nature affected by the resistance variations induced by the PEEP valves. It follows that the flow and concentration values may undergo some variations.

For a more correct flow control at different PEEP values, it is recommended to refer to the present document.

Please note that flowmeter A is used for functioning of the Venturimeter:

+ flow in A = + flow to the patient

Flowmeter B, on the other hand, allows increasing the concentration of oxygen:

+ flow in B = + FI_{O_2} .

Precautions for use of the flowmeters

After using the device, it is essential to turn the flowmeter stopcocks until they are fully closed so that gas delivery is interrupted only after having removed Ventukit Next from the patient.

For use of the flowmeters, always refer to the individual instructions for use contained in the product package.

100%

To deliver 100% FI_{O_2} occlude the air inlet of the Venturi with the blue collar. Connect both A and B lines to an oxygen source and provide a total flow of at least 40 L/min.

Refer to the table attached to the device.

Warnings

It is recommended to use the settings in the table provided by Intersurgical.

These settings assure sufficient flows for proper CPAP therapy with helmet. Higher flow settings would generate unnecessary oxygen consumption and excessive system noise. Should parameters different from those listed in the table be set, it is recommended to check the therapy conditions with the appropriate instruments (oximeters, flowmeters, etc.).

HOW TO ADMINISTER 100% OXYGEN TO THE PATIENT


Fig. 3

the flowmeters A and B.

Fig. 4

If 100% FI_{O_2} needs to be supplied, it is essential to deactivate operation of the Venturimeter.

This is achieved by sliding the blue ring nut (F) onto the air inlet of the Venturimeter (ref. Figs. 3 and 4) The total flow to the patient will be the sum of the flows indicated by

To restore operation of the Venturimeter and thus reduce the FI_{O_2} delivered, slide back the ring nut as shown in Fig. 3 restoring the therapy conditions according to the present document.

Warnings

Always check that the ventilation circuit is perfectly connected.

In the event that the tubes disconnect, patient ventilation will fail. It is recommended to frequently monitor the state of the patient during the therapy.

To assure that the FI_{O_2} delivered to the patient is accurately measured, connect an oximeter to the monitoring port and run random tests.

A. SETTING UP THE DEVICE

Remove the device from the package and model it with one hand in such a way that its cylindrical shape is restored. Open the access port before placing the helmet on the patient. Pre-check the valve. Pull and release the knob checking that the components slide smoothly.

1. Connect the lines A and B to flow as explained above.

It is recommended to set a PEEP of at least 5cm/H₂O as prescribed by the physician. The integrated manometer can be useful to check if there is pressure in the helmet (the indication is however approximate). For the versions not equipped with an integrated manometer an external one can be connected to the pressure port or to an airtight access.

2. Check that the airtight accesses are correctly positioned. Set up the parameters as for the included table, according to physician's prescription.

Activate the ventilation flow before continuing with the next positioning phase in order to reduce the time required to pressurise the helmet.

3. During extended CPAP therapies, the noise inside the helmet might create some discomfort to the patient. To prevent this, Ventukit NEXT is equipped with an integrated silencer that can reduce the noise generated by the Venturimeter by approximately 10 dB. The kit also contains earplugs which provide considerable acoustic comfort for the patient.

B. POSITIONING AND ACTIVATING THE DEVICE

1. Open out the collar (P) using 4 hands in such a way that the patient's head can pass through. To facilitate fitting the helmet, it is suggested that two people open the collar by gripping the flaps in diametrically opposite way and holding their thumbs on the rigid ring of the helmet.

2. The helmet has airtight accesses to insert probes/catheters with diameters between 3.5 and 7 mm.

3. **Access next to the port:** Remove the access from its rigid seat and open it. Position the probe/catheter and refit the access in its rigid seat. Alternatively, remove only the cap and use the central hole. Use this access for NGT if in use.

4. Position the inflatable cushion (O) around the patient's neck and hook the fasteners (M) onto the knobs (I).

5a. Take the rigid annular body (N), open it and position it around the patient's neck following the indications on the body.

5b. Close the rigid annular body on the back of the patient's neck.

6. Using the fastening tabs (L) hook the rigid annular body (N) onto the knobs (I) positioned on the helmet ring.

7. Check that the cushion inflation line (R) runs underneath the rigid annular body (N) and fit it into place in the tube seat (Q) in the lower part of the rigid annular body.

8. Use a pump (S) to inflate the cushion (O) through the inflation line (R) and close the clamp. The cushion serves to improve patient comfort and to hold the helmet in place during the treatment under pressure. Make sure that the elastic collar (P) juts out of the central hole of the cushion (O).

9. Close the access port (C). Pressurise the system by pulling the knob of the anti-asphyxiation valve until the helmet reaches the desired internal pressure. The valve reverts to a safe state to protect the patient when, for any reason, the internal helmet pressure drops to below approximately 2 cmH₂O.

10. For proper functioning of the valve, do not obstruct the valve area and do not lock the action mechanism.

C. USING THE ACCESS PORT (when available)

1. Open the access port to the patient by twisting off the anti-asphyxiation valve anticlockwise.

2. Access the patient's face.

3. Close the access port to the patient by twisting on the anti-asphyxiation valve clockwise.

4. Pull the valve knob using two fingers until the helmet is pressurised.

D. HOW TO REMOVE THE HELMET

1. Remove any probes/catheters used during the therapy from the airtight accesses (D). Open the clamp of the inflation line (R) and let the cushion deflate holding one hand on the top of the helmet.

2. Open the access port.

3. Remove the cushion inflation line (R) from its seat (Q). Undo the fastening tabs of the rigid annular body (N) thus releasing the helmet. Undo the cushion fasteners (M) from the front knobs on the helmet.

4. Move the deflated cushion to the back of the patient's neck. Remove the helmet by opening out the collar with 4 hands. Also remove the rigid annular body if it is no longer needed.

5. Turn off the ventilation flow.

WARNINGS / PRECAUTIONS

 Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

 Caution: this product contains phthalate. There is evidence suggesting that

phthalate exposure from medical procedures may harm children, pregnant women or breastfeeding women. Use of this product does not represent a high exposure risk to phthalate, however, as a precautionary measure its use should be limited to essential procedures.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.

2. Monitor the clinical parameters of the patient. Warning: models without anti-asphyxiation valve must only be used if monitoring of the clinical parameters is guaranteed with adequate systems equipped with alarms.

3. If the patient's conditions DETERIORATES or there is no improvement within the scheduled times, it is recommended to evaluate alternative ventilation Techniques.

4. The efficacy of the therapy is considerably affected by the PEEP value set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient.

Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.

5. Caution: the device contains metal; do not use in departments where the presence of metal is a source of risk to the safety of the patient or third parties.

6. During therapy, the patient may complain of annoying noise caused by the high flows utilised. This is due to the velocity and turbulence of the gas. The annoying noise can easily be reduced by fitting one or two combined filters for breathing circuits on the hood connectors (normally found in hospital departments and not provided in this kit) or having the patient wear sound-absorbing ear plugs (provided in this kit).

7. When used with humidified ventilation gases, check that any condensate forming in the helmet does not cause the patient discomfort.

8. The patient must REMOVE all accessories such as earrings, hair clips, combs and any other metal objects.

9. The physician is responsible for deciding which ventilation therapy is most suited to the patient's pathology. The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. The neck circumferences (expressed in centimetres) and the helmet size is indicated on each device. Suitable choices will ENSURE A GOOD SEAL around the patient's neck.

10. Use ONLY UN BROKEN skin.

11. The device can be used for maximum 7 days. Depending on the clinical picture of the patient, the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.

12. Non-sterile. Do not sterilize.

13. Do not use on another patient.

14. Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°/+50°C).

15. To facilitate modeling/shaping the helmet, it is recommended to keep it at room temperature (about 20- 25°C) for a few hours before use.

16. Dispose of the materials immediately after use CONFORMING TO the current laws and regulations.

17. Do not use pressurised gas sources to inflate the integrated cuff.

18. Activating the anti-asphyxiation valve in the event of ventilation interruption allows air exchange with the outside limiting CO₂ rebreathing. The anti-asphyxiation valve does not substitute the ventilation support in the event of accidental ventilation interruption. **The anti-asphyxiation valve does not exempt the nursing staff from implementing appropriate monitoring and supervision.**

19. Caution: Use of the helmet does not exempt operators from the obligation to use adequate Personal Protection Devices in accordance with the hospital procedures.

20. Non invasive ventilation shall be interrupted if there is one of the following conditions:

- worsening state of consciousness and respiratory distress
- loss of airway protection
- unchanged PaCO₂ (during two subsequent EGA performed after max 1 hour)
- persistence of severe hypoxemia
- severe and uncontrollable hemodynamic instability
- patient / ventilator synchronization problems
- uncontrollable secretions
- interface intolerance.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy