VENTUKIT Up

- it** CASCO PER N.I.V. (CPAP) CON SISTEMA DI MISCELAZIONE ARIA/OSSIGENO
- en** HELMET FOR N.I.V. (CPAP) WITH AIR/OXYGEN BLENDING SYSTEM
- de** HELM FÜR N.I.V. (CPAP) MIT SYSTEM ZUR LUFT-SAUERSTOFF-MISCHUNG
- sr** KACIGA ZA N.I.V. (CPAP) SA SISTEMOM ZA MEŠANJE VAZDUHA / KISEONIKA
- tr** NIV (CPAP) İçin Hava/oksijen karışım sistemli Başlık

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INSTRUCTIONS FOR USE

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
UPUTSTVO ZA UPOTREBU

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

E

F1

F2

F3


F4

F6


B1


B2/B3

B3a


B4

B5


B6


B7


B7a


B7b


B8


B9


B11


B12


C1


C2

C3

C4

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KIT CONTENTS

- Non-invasive ventilation hood (CPAP) with integrated air/oxygen mixing system.
- Oxygen feed tubes
- PEEP valve and pressure gauge
- Anti-noise earplugs
- Instructions for use.

Note: the flowmeters are not included in the kit.

Intended use

Ventukit UP is a device intended for non-invasive CPAP ventilation (Continuous Positive Airway Pressure) in a hospital environment when an FiO_2 variable between 40% and 100% needs to be supplied.

Indications for use

Suitable for treatment of patients with acute hypoxemic respiratory failure. For other pathologies the treatment may only be performed under careful monitoring by expert operators.

Contraindications

- Coma
- Uncooperative patient
- Cardiac arrest
- Haemodynamic instability
- Recent oesophageal and gastro-surgical operations
- Heavy bleeding of the upper digestive tract
- Obstruction of the upper airways
- Pneumothorax

Restrictions of use

1. The device must be used by qualified and trained medical/nursing staff.
2. To be used for administration of air and oxygen.
3. For pure oxygen treatment, connect both tubes to oxygen sources and feed the device with O_2 flows of not less than 40 l/min deriving from the sum of the two lines.
4. If used properly, the device can be used continuously for a maximum of 7 days after which it must be replaced.
5. Monitor the clinical parameters of the patient.
6. To be used solely with oxygen flowmeters conform to the standards and capable of delivering the flows indicated in the table. If fed with an O_2 source (only line A), Ventukit UP can reach an FiO_2 of 40%. To achieve a higher FiO_2 use two O_2 sources at the same time (line A and line B). The use of flow meters with an output of less than the maximum required in the table on line B, will reduce the performance of the device.
7. Use flowmeters with a 7mm diameter hose tail.
- 8.

PATIENTS	RECOMMENDED SIZE
ADULTS - weighing > 30 Kg	ADULT XS-S-M-L-XL-XXL, based on the neck circumference indicated on the label

The medical/nursing staff is responsible for choosing the

device most suitable for the patient in terms of configuration and size. An appropriate choice will ensure good hermetic seal around the patient's neck. The neck circumferences are approximate; the physician is advised to measure the neck circumference of the patient to be treated.

It is recommended to use the included measuring tape to choose the best size according to the patient's neck circumference.

Where and by whom can Ventukit UP be used

Ventukit UP may only be used by qualified medical/nursing staff in a hospital environment. The therapy is started simply by connecting Ventukit UP to the oxygen source available in hospital wards (wall-mounted distribution systems or cylinders), for example: emergency medicine, pneumology, haematology, oncology, post-surgery intensive care, ambulance, emergency room, etc.

Complications and side effects

Episodes of armpit pain may occur in patients that have a sensitive skin and during extended therapies at high PEEP pressures, which can be avoided by applying thin strips of hydrocolloid or by alternating the fastening systems such as chest straps or elastic bands for fastening to the bed. Some patients may suffer episodes of claustrophobia. The helmet is equipped with a quick-access opening so that the operator can easily access the patient's face. Improper use of the device resulting in flows of less than 40 l/min may cause an increase of CO_2 inside the helmet.

Ventukit UP, if used in compliance with these instructions, guarantees a total flow to the patient of more than 40 l/min thus obtaining proper ventilation gas exchange.

During extended CPAP therapies, the noise inside the helmet might create some discomfort to the patient. To prevent this, Ventukit UP is equipped with an integrated silencer that can reduce the noise generated by the Venturimeter by approximately 10 dB. The kit also contains earplugs which provide considerable acoustic comfort for the patient.

How to use the VENTUKIT UP

The therapy is started simply by connecting Ventukit UP to the available oxygen source.

The Ventukit UP tubes A and B then need to be properly connected to the respective flowmeters (Ref. Figs. 1 and 2) connected to the oxygen source in the ward appropriately pressurised following the instructions provided by the manufacturer of the flowmeters.

Ventukit UP is provided with a performance table of the system's operating conditions.

Example: to administer 50% FiO_2 + a PEEP of 7.5 cm H_2O set the flowmeter A to 8 l/min and the flowmeter B to 9 l/min, thus obtaining a flow of about 47 l/min to the patient.

The performances of the Venturimeters are by their nature affected by the resistance variations induced by the PEEP valves. It follows that the flow and concentration values may undergo some variations. For a more correct flow control

at different PEEP values, it is recommended to refer to the present document.

Please note that flowmeter A is used for functioning of the Venturimeter: **+ flow in A = + flow to the patient**

Flowmeter B, on the other hand, allows increasing the concentration of oxygen: **+ flow in B = + FiO₂**

Fig. 1


Fig. 2

Precautions for use of the flowmeters

After using the device, it is essential to turn the flowmeter stopcocks until they are fully closed so that gas delivery is interrupted only after having removed Ventukit UP from the patient.

For use of the flowmeters, always refer to the individual instructions for use contained in the product package.

100%	To deliver 100% FiO ₂ occlude the air inlet of the Venturi with the blue collar. Connect both A and B lines to an oxygen source and provide a total flow of at least 40 L/min.
-------------	---

Refer to the table attached to the device.

Warnings

It is recommended to use the settings in the table provided by Intersurgical.

These settings assure sufficient flows for proper CPAP therapy with helmet. Higher flow settings would generate unnecessary oxygen consumption and excessive system noise. Should parameters different from those listed in the table be set, it is recommended to check the therapy conditions with the appropriate instruments (oximeters, flowmeters, etc.).

HOW TO ADMINISTER 100% OXYGEN TO THE PATIENT

Fig. 3

If 100% FiO₂ needs to be supplied, it is essential to deactivate operation of the Venturimeter.

This is achieved by sliding the blue ring nut "G" onto the air inlet of the Venturimeter (ref. Figs. 3 and 4).

The total flow to the patient will be the sum of the flows indicated by the flowmeters A and B.

Fig. 4

To restore operation of the Venturimeter and thus reduce the

FiO₂ delivered, slide back the ring nut as shown in Fig. 3 restoring the therapy conditions according to the present document.

Warnings

Always check that the ventilation circuit is perfectly connected. In the event that the tubes disconnect, patient ventilation will fail. It is recommended to frequently monitor the state of the patient during the therapy.

To assure that the FiO₂ delivered to the patient is accurately measured, connect an oximeter to the monitoring port and run random tests.

A. SETTING UP THE DEVICE

Remove the device from the package and model it with one hand in such a way that its cylindrical shape is restored.

Open the access port (if present) before positioning the helmet on the patient. Pre-check the valve. Pull and release the knob checking that the components slide smoothly.

1. Separate the two parts of the device (upper and lower rings).
2. Arrange the fastening systems securing them on the rear of the helmet.

B. POSITIONING AND ACTIVATING THE DEVICE

1. Positioning of the lower ring.

Open out the collar using 4 hands, in such a way that the patient's head can pass through.

To facilitate fitting the helmet, it is suggested that two people open the collar by gripping the flaps in a diametrically opposite way and holding their thumbs on the rigid ring of the helmet.

2. Complete positioning of the fastening system on the front of the Helmet.
3. The helmet has airtight accesses to insert probes/catheters with diameters between 3.5 and 7 mm. Pass the end connectors of any probes or catheters through the appropriate accesses on the collar of the helmet from the inside to the outside of the membrane.
 - Remove the access from its rigid seat and open it
 - Position the probe/catheter and refit the access in its rigid seat
 - Alternatively, remove only the cap and use the central hole
4. Connect the lines A and B to flow as explained above.

Connect the PEEP valve to the outlet connector.

It is recommended to set a PEEP of at least 5cm/H₂O as prescribed by the physician. The integrated manometer can be useful to check if there is pressure in the helmet (the indication is however approximate).

For the versions not equipped with an integrated manometer an external manometer can be connected to the pressure port or to an airtight access.

5. Check that the connection is secure before proceeding with administration of the therapy. Activate the ventilation flow before continuing with the next positioning phase in order to reduce the time required to pressurise the

helmet. Set up the parameters as for the included table, according to physician's prescription.

6. During extended CPAP therapies, the noise inside the helmet might create some discomfort to the patient. To prevent this, Ventukit Up is equipped with an integrated silencer that can reduce the noise generated by the Venturimeter by approximately 10 dB. The kit also contains earplugs which provide considerable acoustic comfort for the patient.
7. **Positioning of the upper ring.** Attach the upper part of the helmet to the lower ring, previously placed on the patient, and lock into place using the levers.
9. Pressurise the system by pulling the knob of the anti-asphyxiation valve until the helmet reaches the desired internal pressure. The valve will go into patient safety and protection mode when, for any reason, the internal helmet pressure drops below 2 cm H₂O.
10. For proper functioning of the valve, do not obstruct the valve area and do not lock the action mechanism.
11. Adjust the length of the straps so that the rigid ring is about 1 centimetre from the patient's shoulders.
12. Some models have an internal inflatable cuff.

C. USING THE ACCESS PORT (WHEN AVAILABLE)

1. Open the access port to the patient by turning the anti-asphyxiation valve anticlockwise.
2. Access the patient's face.
3. Close the access port to the patient by turning the anti-asphyxiation valve clockwise.
4. Pull the valve knob using two fingers until the helmet is pressurised.

D. TO PERFORM STANDARD NURSING PROCEDURES OR FOR EMERGENCY ACCESS TO THE PATIENT

1. Open the access port to the patient (if present).
2. Separate the upper part of the helmet from the lower ring using the levers located on the rigid rings. Perform procedures on the patient.
3. Reposition the upper part of the helmet onto the lower ring remaining on the patient, and lock into place using the levers.
4. Close the access port (if present).
5. Pressurise the system by pulling the knob of the anti-asphyxiation valve until the helmet reaches the desired internal pressure.

E. SYSTEM TO IMPROVE COMFORT

Alternative fastening systems

During extended therapy and with high operating pressures the patient may complain of armpit pain. Intersurgical can supply different accessories to give the patient relief:

- An abdominal/chest belt to be fixed
- An elastic strap with slots to fix the helmet to the bed.

To alleviate pain, hydrocolloid dressings may be applied between the patient's skin and the axillary protection.

F. HOW TO REMOVE THE HELMET

1. Deflate the internal cuff (where present) by opening the clamp. Open the access port to the patient (if present).
2. Separate the upper part of the helmet from the lower ring using the levers located on the rigid rings.
3. Turn off the ventilation flow.
4. Remove any probes/catheters used during the therapy from the air tight accesses.
5. Detach the fastening systems.
6. Remove the lower part of the helmet. Open out the collar using 4 hands, in such a way that the patient's head can pass through.

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-  Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

 Caution: this product contains phthalate. There is evidence suggesting that phthalate exposure from medical procedures may harm children, pregnant women or breastfeeding women. Use of this product does not represent a high exposure risk to phthalate, however, as a precautionary measure its use should be limited to essential procedures.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
2. Monitor the clinical parameters of the patient.
3. If the patient's conditions DETERIORATES or there is no improvement within the scheduled times, it is recommended to evaluate alternative ventilation techniques.
4. The efficacy of the therapy is considerably affected by the PEEP values set. Therefore, the medical staff are advised to carefully evaluate the pressure level most suited to the clinical condition of the patient.
Using too low pressure may not be sufficient for alveolar recruitment. Using too high pressure may be a source of alveolar over distension.
5. Caution: the device contains metal; do not use in departments where the presence of metal is a source of risk to the safety of the patient or third parties.
6. The annoying noise can easily be reduced by fitting one or two combined filters for breathing circuits on the hood connectors (normally found in hospital departments and not provided in this kit) or having the patient wear sound-absorbing ear plugs (provided in this kit).
7. When used with humidified ventilation gases, check that any condensate forming in the helmet does not cause the patient discomfort.
8. The patient must REMOVE all accessories such as earrings, hair clips, combs and any other metal objects.
9. When fastening the device, make sure that no parts, such as hair or fabrics, become trapped between the upper and

- lower rings.
10. The physician is responsible for deciding which ventilation therapy is most suited to the patient's pathology. The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. The neck circumferences (expressed in centimetres) and the helmet size is indicated on each device. Suitable choices will ENSURE A GOOD SEAL around the patient's neck.
 11. Use ONLY ON UNBROKEN skin.
 12. The device can be used for maximum 7 days.
Depending on the clinical picture of the patient, the medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
 13. Non-sterile. Do not sterilize.
 14. Do not use on another patient.
 15. Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (20°/+50°C).
 16. To facilitate modelling/shaping the helmet, it is recommended to keep it at room temperature (about 20-25°C) for a few hours before use.
 17. Dispose of the materials immediately after use CONFORMING TO the current laws and regulations.
 18. Activating the anti-asphyxiation valve in the event of ventilation interruption allows air exchange with the outside limiting CO₂ rebreathing. The anti-asphyxiation valve does not substitute the ventilation support in the event of accidental ventilation interruption. **The anti-asphyxiation valve does not exempt the nursing staff from implementing appropriate monitoring and supervision.**
 19. Caution: Use of the helmet does not exempt operators from the obligation to use adequate Personal Protection Devices in accordance with the hospital procedures.
 20. The cuff inflation pump could not be included in the kit. Do not pressurised gas sources to inflate the integrated cuff.
 21. Non invasive ventilation shall be interrupted if there is one of the following conditions:
 - worsening state of consciousness and respiratory distress
 - loss of airway protection
 - unchanged PaCO₂ (during two subsequent EGA performed after max 1 hour)
 - persistence of severe hypoxemia
 - severe and uncontrollable hemodynamic instability
 - patient / ventilator synchronization problems
 - uncontrollable secretions
 - interface intolerance

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Intersurgical SpA - Via Morandi, 12 - 41037 Mirandola (MO) Italy