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en **Instructions for use - Reusable mask for non-invasive ventilation (NIV) –
with reusable connectors**

INTENDED USE

Reusable mask for CPAP or BILEVEL ventilatory support.

INDICATIONS FOR USE

Indicated for the treatment of respiratory failure using non-invasive CPAP or BILEVEL ventilation methods.
The mask can be used with more than one patient, both adult and pediatric ones.

CONTRAINDICATIONS

- The device is not to be used on uncooperative or unconscious patients.
- The device is not suitable for use on patients that are unable to remove the mask on their own, take drugs that may cause vomiting, have night-time vomiting, hiatus hernia or functional damage to the cardiac sphincter.
- This mask is not suitable for use on patients that have the following conditions: glaucoma, recent eye surgery or eye dryness.
- This mask cannot be used for uncooperative patients and unable to remove the mask on their own.
- This mask is not suitable for the administration of anesthetic gases and drugs.

COMPLICATIONS/SIDE EFFECTS

Unusual chest discomfort, wheezy breathing, gastric dilatation, belching, severe headache, eye dryness or pain, eye infections, blurred vision during use of the positive pressure device must be reported to the attending physician.
Claustrophobia.
CO2 retention (rebreathing).
Facial irritations and traumas.
Skin irritations
Dry eyes and airways
pneumothorax

WARNINGS AND PRECAUTIONS

In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
Monitor the vital parameters of the patient during the treatment.
If noticing that the patient's conditions worsen or that there is no improvement within the scheduled times, it is recommended to evaluate using an alternative ventilation technique.
When used with humidified ventilation gases, check that any condensate forming in the mask does not cause the patient discomfort.
The medical/nursing staff is responsible for choosing the device most suitable for the patient in terms of configuration and size. The appropriate choice will guarantee good seal.
Use on undamaged skin.
Expiry: 5 years if the packaging is intact and if stored in normal storage conditions.

Non-invasive ventilation shall be interrupted if there is one of the following conditions:

- worsening state of consciousness and respiratory distress
- loss of airway protection
- unchanged PaCO2 (during two subsequent EGA performed after max 1 hour)
- persistence of severe hypoxemia
- severe and uncontrollable hemodynamic instability
- patient / ventilator synchronization problems
- uncontrollable secretions
- interface intolerance

For the use of connectors without asphyxiation valve:

Use only with fans with adequate alarms and safety systems that detect ventilation failures. This mask does not include anti-asphyxiation valves that would allow the patient to breathe even in the event of a ventilator failure.

Mask with standard connector (BLUE) or bronchoscopy access connector (GREEN)	Mask with antisphyxia valve connector (ORANGE)
<p>This mask does not have an integrated air outlet port and the circuit hence requires a separate expiration device. Use the expiration port recommended by the device manufacturer.</p> <ul style="list-style-type: none"> • This device is not suitable for ventilatory support of the vital functions. • To reduce the risk of vomiting in sleep to a minimum, the patient should not eat or drink three (3) hours before use. • Check that the mask and the cranium support strap are the right size for the patient's face. If the patient is very uncomfortable, rearrange the mask and the cranium support strap. If the mask is still uncomfortable, use another mask. 	<p>This mask has an integrated air outlet port which provides a continuous patient circuit leakage path. It is recommended that the circuit have a separate expiration device. It may be necessary to adjust the pressure level and the oxygen therapy to compensate for the additional leakage created by the expiration device.</p> <ul style="list-style-type: none"> • This device is not suitable for ventilatory support of the vital functions. • To reduce the risk of vomiting in sleep to a minimum, the patient should not eat or drink three (3) hours before use. • Check that the mask and the cranium support strap are the right size for the patient's face. If the patient is very uncomfortable, rearrange the mask and the cranium support strap. If the mask is still uncomfortable, use another mask. • This CPAP mask may exclusively be used with CPAP systems according to the instructions of the attending physician or the pneumologist. Use the mask only when the CPAP system is on and functioning properly. The breather and other holes in the mask must never be blocked. <p><i>Why this warning? CPAP systems are designed to be used with special masks equipped with connectors with breather holes that allow a continuous air flow outside the mask. When the CPAP device is on and working, the fresh air delivered by the CPAP device pushes the exhaled air to the outside through the mask expiration port. However, when the CPAP device is off, insufficient fresh air is delivered and consequently the exhaled air could be re-inhaled. In some circumstances, if closed-circuit breathing is prolonged for several minutes, it can cause suffocation. These warnings refer to the majority of CPAP systems.</i></p>

LIMITS OF USE

Use exclusively in hospital or healthcare environments.
The device must be used by qualified and trained medical/nursing staff.
Use with monitoring lines suitable to be connected to pressure taps with OD 5 mm.
SNG cap (if present) can be used to introduce nasogastric tubes with diameter ≥ 12 Ch.

WARNINGS RELATED TO REUSE

The reconditioning instructions contained in this package insert must be followed before each subsequent reuse of the device.

Do not use gamma rays or ultrasonic cleaning equipment. Do not sterilize with dry or wet steam sterilizers. Do not sterilize or disinfect with components or derivatives of phenol and alkylamine (glucoprotamine).

Do not clean or rub with solvent or chloride based solutions. The necks must be removed from the masks before disinfection and / or sterilization.

The connectors must be removed before disinfection and / or sterilization.

Do not use hard brushes or other materials that could damage the silicone parts of the device. Avoid contact with hard or sharp objects.

After cleaning, disinfection and sterilization, visually inspect the product. Slight yellowing / color changes are considered normal. Do not reuse in case of evident cracks or detachments of the components.

In the event that the products, despite the cleaning procedure, are not cleaned properly, it is not possible to proceed with their disinfection or sterilization, but the aforementioned devices must be disposed of as per the dedicated paragraph.

After disinfection and sterilization, keep the mask properly so as not to pollute it.

Follow the technical data sheets of the disinfectants for operating procedures and safe handling of the same.

LIMITS OF REUSE

Device sold in non-sterile configuration. Reusable and reconditionable up to a maximum of 10 times.

The device can be reprocessed by high level disinfection (CIDEX OPA for mask and fittings, ANIOXIDE 1000 only for the mask) and / or sterilization (ETO and / or Hydrogen peroxide plasma gas).

INSTRUCTIONS FOR USE

INSPECTION AND FUNCTIONAL TEST BEFORE USE

- Carefully read the instructions.
- Before each use and before the first use, check that there are no damages or deteriorations that could affect the correct functioning of the device. Check the connectors for any deformations, breakages, leaks and obstructions.
- Verify the firm connection between the elbow and the mask shell.
- Thoroughly clean the patient's face each time before use.
- Choose the correct mask size: the mask gasket must be positioned above the eyebrow and directly under the lower lip (refer to figure 1).
- Place the mask on the patient's face adjust the straps by pulling on the fixing neck (refer to figure 2).

Figure 1


Figure 2

For models equipped with anti-asphyxiation valve:

- Check functioning of the safety valve on the anti-asphyxiation elbow.
- When the unit is off, the safety valve membrane must be open to allow the patient to breathe in ambient air.
- Turn on the CPAP or BILEVEL unit.
- The hinge valve must close so that the patient receives air from the unit.
- If the valve does not work as indicated, reposition the mask.
- Do not obstruct the safety valve opening.
- Always check the valve before using it to identify any malfunctions or obstructions caused by patient secretions.

Use of the interchangeable elbow accessory (ref. codes in the table)

Code	Description	Colour
MASKR-ACC9599	Standard angled connector, 22F connection - reusable	Blue
MASKR-ACC9593	Angled connector with anti-asphyxiation and outlet valve, 22F connection - reusable	Orange
MASKR-ACC9597	Angled connector with bronchoscopy access, 22F connection - reusable	Green

Remove the mask from the patient.

Remove the existing elbow from the mask holding the adapter ring and gently pulling and turning the elbow.

Gently push the desired compatible interchangeable elbow into the adapter ring of the mask until it perfectly fits into place. Do not force it or use tools.

Use of the standard interchangeable elbow accessory for bronchoscopy

The standard interchangeable bronchoscopy connector is an accessory for a full facemask that allows delivering the CPAP and BILEVEL therapy to the patient in the course of a bronchoscopy procedure. The use of the elbow with this mask is an option with CPAP and BILEVEL systems for the treatment of obstructive sleep apnoea syndrome (OSAS) and/or ventilatory support.

NOTE

The standard interchangeable bronchoscopy connector is to be used temporarily only during bronchoscopy procedures.

CONTRAINDICATIONS

- The device is not to be used on uncooperative or unconscious patients.
- The device is not suitable for use on patients that are unable to remove the mask on their own and who take drugs that may cause vomiting.

WARNINGS

- This connector may be used for treatments with a nebulizer exclusively in the presence of trained medical staff.
- Use only with a ventilator equipped with appropriate alarms and safety systems against any ventilator malfunctioning. This mask does not include an anti-asphyxiation valve, which would allow the patient to breathe even in the event of ventilator malfunctioning.
- During use of the bronchoscopy connector, the mask does not have an integrated air outlet port and the circuit hence requires a separate expiration device. Use the expiration port recommended by the device manufacturer.
- This device has been designed to be used exclusively with a mask compatible with the interchangeable elbow of the full facemask.

INSTRUCTIONS FOR RECONDITIONING

POINT OF USE

Please refer to cleaning instructions below.

CONTENT AND TRANSPORTATION

It is recommended to recondition all products, as soon as reasonably possible, after use and before subsequent reuse.

PREPARATION BEFORE CLEANING

Before proceeding with cleaning and disinfection or sterilization, separate the interchangeable connector and nugal piece from the mask.

The device for fastening the mask on the patient's head may not be treated with the disinfection and sterilization methods listed, in fact the nugal piece is single-patient accessory and must always be replaced before reusing the mask.

The SNG cap (if present) is a single patient accessory and must always be replaced before the mask can be reused.

A slight smell after disinfection is normal.

Cleaning is an obligatory step before any subsequent reconditioning. Before using, disinfecting and sterilizing, clean the mask and the fittings according to the indications below.

Manual cleaning:

Wash your hands with a mild detergent such as a neutral pH liquid soap or trisodium phosphate. Rinse thoroughly with clean water, and dry with a mild absorbent material that does not release fibers. Disassemble the mask by removing the neck, elbow fitting and SNG cap (if present).

Open the pressure plugs of the mask.

Immerse the device for 5 minutes in a commercially available anionic detergent, such as DIALZIMA ULTRA diluted in water (0.15%). Rinse for 10 seconds under running water.

Do not use hard brushes or other materials that could damage the products. Do not use a washing machine. Dry at room temperature. Do not dry on a hot surface. Dry in the air avoiding direct exposure to sunlight.

WARNING

After cleaning, disinfection or sterilization can be carried out before reuse.

Before each subsequent disinfection or sterilization cycle, always proceed with the preparation operations before cleaning and cleaning itself listed above.

High level manual disinfection

Treatment	Exposure time	Active ingredient	Solution	Rinsing	Maximum disinfection cycles
CIDEX OPA	5 min. immersion	0.55% ortho-phthalaldehyde	/	10 min. immersion in water (3 times)	10
ANIOXIDE 1000 (just for mask and not for connectors)	5 min. immersion	0.15% peracetic acid	Activation as per the manufacturer's technical data sheet	1 min. immersion in water (once)	10

Sterilization:

Intersurgical SpA verified the effectiveness of sterilization on medical devices named MASKR. The parameters below are to be considered recommended. Any different cycles must be validated by the customer / user.

After sterilization, follow the degassing times (for ETO) and finally allow the device to dry at room temperature.

Hydrogen peroxide gas plasma

Treatment	Exposure time	Maximum number of cycles
STERRAD	See the machine manual	10

ETO

Treatment	Exposure time	Maximum number of cycles
Ethylene Oxide	18 hours +6 degassing hours	10

STORAGE

The device must be stored at a temperature between -20 ° C - +50 ° C. Keep away from direct light, from sources of radiation, fluorescent light, electric motors and diathermy devices. The degradation of the rubber increases in the presence of oxygen, oxidizing agents, ozone, ether, phenol, cresol, terpene, hypochlorite, formaldehyde, ketone, inorganic acids and hydrocarbons. Avoid contact with metals, oils, solvents, greases, and strong detergents. If rubber products break, crack, lose elasticity, become rigid, or lose their color on contact, they are no longer safe for further use and must be replaced. The instructions provided are validated by the COMPETENT qualified manufacturer in the manufacture of a reusable device. The reconditioner is responsible for the fact that the reconditioning, carried out effectively with equipment, materials and personnel in a specific structure, achieves the desired result. This normally requires validation and continuous monitoring of the process. Likewise, any departure by the reconditioner from the instructions provided must be properly assessed in terms of efficiency and potential adverse consequences.

OPERATING SPECIFICATIONS

+5 - +40°C with 15-95% humidity

DEVICE DISPOSAL

Dispose of the materials as potentially hazardous material in compliance with the legislation in force as regards hospital environments.

ACCESSORIES

Interchangeable elbow reusable connectors and nucal piece (ref. codes in the table)

Code	Description	Colour
MASKR-ACC9599	Standard angled connector, 22F connection - reusable	Blue
MASKR-ACC9593	Angled connector with anti-asphyxiation and outlet valve, 22F connection - reusable	Orange
MASKR-ACC9597	Angled connector with bronchoscopy access, 22F connection - reusable	Green
313-8344	Replaceable nucal piece for MaxShield Mask – Large / Extra Large Mask size	/
313-8347	Replaceable nucal piece for MaxShield Mask – Small size	/

AS THE DEVICE IS REUSABLE, KEEP THESE INSTRUCTIONS FOR USE

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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