
- it** Accessorio a boccaglio per casco
- en** Mouthpiece accessory for helmet
- fr** Accessoire à embout per casque
- de** Mundstück-Zubehör für Helm
- es** Accesorio boquilla para casco
- pt** Bocal para capacete
- tr** Başlık için ağızlık aksesuarı

- nl** Accessoire met mondstuk voor kap
- sv** Anslutningskopplingstillbehör för hjälm
- el** Επιπτόμο για κάσκα
- sl** Priključek ustnika za skafander
- ru** Принадлежность для шлема - мундштук
- sk** Príslušenstvo s náustkom pre helmu
- sr** Usni dodatak za kacigu

**Istruzioni per l'uso / Instructions for use / Mode d'emploi / Gebrauchsanleitung / Instrucciones de uso / Instruções de utilização / Gebruiksaanwijzing / Bruksanvisning / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / Navodilo za uporabo/ ИНСТРУКЦИЯ К ИСПОЛЬЗОВАНИЮ / Návod na použitie / Kullanma Talimatları / Uputstvo Za Upotrebu**

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**en** INDICATIONS FOR USE

Mouthpiece for aerosol therapy to use in combination with Intersurgical helmets with patient access port.

**RESTRICTIONS OF USE**

1. Connectable to therapeutic respiratory device fitted with 22mm Female connection (22 mm ID).
2. The device must be used with Intersurgical helmets for CPAP and NIV therapy with patient access port.
3. The device must be used by qualified and trained medical/nursing staff.

**INSTRUCTIONS FOR USE**

- A) Remove the device from the package;
- B) assemble the elbow connector 4 and the mouthpiece 5 on the screwed cap 3 as shown in the drawing;
- C) remove the patient access port fitted with bi-directional anti-asphyxiation valve (2) from the helmet (1) by turning counterclockwise;
- D) connect the therapeutic respiratory device to the elbow adapter 4;
- E) verify the security of connections;
- F) the device is ready to use. To deliver therapy, follow the instructions for use provided by the therapeutic device manufacturer;
- G) At the end of the aerosol therapy it is necessary to remove the mouthpiece accessory and reattach the patient access port fitted with the bi-directional anti-asphyxiation valve (2) to the helmet (1) by turning clockwise.

**WARNINGS / PRECAUTIONS**

Single use device. Reuse is not allowed since it could cause cross contamination in patients. Moreover, materials which make part of the device could be damaged after cleaning/disinfection and therefore they could not guarantee required performances, therefore creating a safety risk for patients.

1. In the event of an evident change in performance or inadequate performance of the device, it is recommended to replace it.
2. Monitor the clinical parameters of the patient during use. Since the bi-directional anti-asphyxiation valve is not in place during treatment, close monitoring of the patient is required.
3. The medical/nursing staffs are responsible for defining the need for more frequent replacement of the device.
4. Do not attach ventilator tubing to the elbow connector 4.
5. Non-sterile. Do not sterilize.
6. Do not use on another patient.
7. Expiry: 5 years provided that the packaging is undamaged and if stored in normal storage conditions (-20°C/+50°C).
8. Dispose of the materials immediately after use conforming to the current laws and regulations.
9. If any visible deterioration of the mouthpiece is apparent (cracking, crazing, broken, tears etc), it should be discarded and replaced.
10. Using a mouthpiece may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult your physician or dentist.

**Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

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