




Intersurgical Venturi valves – Minimum flow rates




Product Code: 1400024
Venturi valve 24% oxygen
Blue
2 L/min




Product Code: 1400028
Venturi valve 28% oxygen
White
4 L/min




Product Code: 1400031
Venturi valve 31% oxygen
Orange
6 L/min




Product Code: 1400035
Venturi valve 35% oxygen
Yellow
8 L/min



Product Code: 1400040
Venturi valve 40% oxygen
Red
10 L/min



Product Code: 1400050
Venturi valve 50% oxygen
Pink
15 L/min



Product Code: 1400060
Venturi valve 60% oxygen
Green
15 L/min