



AccuPAP™


Positive Airway Pressure Therapy System



Oxygen and Aerosol Therapy ▪ Pulmodyne®

AccuPAP™

The AccuPAP system is designed to deliver positive airway pressure indicated for the prevention or reversal of atelectasis in adult patients. The AccuPAP includes adjustable PEEP, this allows the healthcare professional to deliver pressure levels according to each patient's requirements. The PEEP can be adjusted without requiring a change to be made to the input gas flow.



Code	Description	Box Qty.
313-6000	AccuPAP with mouthpiece and oxygen tube, 2.1 m	10
313-6001	AccuPAP with mouthpiece, disposable pressure manometer and oxygen tube, 2.1 m	10

[Make an enquiry](#)

[▶ Watch the video](#)

AccuPAP and DPM are manufactured by Pulmodyne Inc, USA and distributed by Intersurgical Ltd

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS10.40 • Issue 7 07.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan
 Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA