



InterTube™ tracheal tubes


A comprehensive range of ET tubes for adult and paediatric use



Airway Management • Airway Devices














InterTube™ tracheal tubes

Designed for use in anaesthesia, A&E and emergency medicine, the InterTube range offers three varieties of tracheal tube: cuffed, reinforced and uncuffed, in a comprehensive range of sizes from paediatric to adult. Supplied sterile, individually packed, latex free and for single use.


















Cuffed

Code	Description	Size (Inner Diameter)	Box Qty.
9040100	InterTube tracheal tube, cuffed	10.0 mm	10 
9040095 	InterTube tracheal tube, cuffed	9.5 mm	10 
9040090	InterTube tracheal tube, cuffed	9.0 mm	10 
9040085	InterTube tracheal tube, cuffed	8.5 mm	10 
9040080	InterTube tracheal tube, cuffed	8.0 mm	10 
9040075	InterTube tracheal tube, cuffed	7.5 mm	10 
9040070	InterTube tracheal tube, cuffed	7.0 mm	10 
9040065	InterTube tracheal tube, cuffed	6.5 mm	10 
9040060	InterTube tracheal tube, cuffed	6.0 mm	10 
9040055	InterTube tracheal tube, cuffed	5.5 mm	10 
9040050	InterTube tracheal tube, cuffed	5.0 mm	10 
9040045	InterTube tracheal tube, cuffed	4.5 mm	10 














Wire-reinforced

Code	Description	Size (Inner Diameter)	Box Qty.
9060100	InterTube tracheal tube, wire-reinforced cuffed	10.0 mm	10 
9060095	InterTube tracheal tube, wire-reinforced cuffed	9.5 mm	10 
9060090	InterTube tracheal tube, wire-reinforced cuffed	9.0 mm	10 
9060085 	InterTube tracheal tube, wire-reinforced cuffed	8.5 mm	10 
9060080	InterTube tracheal tube, wire-reinforced cuffed	8.0 mm	10 
9060075	InterTube tracheal tube, wire-reinforced cuffed	7.5 mm	10 
9060070	InterTube tracheal tube, wire-reinforced cuffed	7.0 mm	10 
9060065	InterTube tracheal tube, wire-reinforced cuffed	6.5 mm	10 
9060060	InterTube tracheal tube, wire-reinforced cuffed	6.0 mm	10 
9060055	InterTube tracheal tube, wire-reinforced cuffed	5.5 mm	10 
9060050	InterTube tracheal tube, wire-reinforced cuffed	5.0 mm	10 
9060045	InterTube tracheal tube, wire-reinforced cuffed	4.5 mm	10 








Uncuffed

Code	Description	Size (Inner Diameter)	Box Qty.
9050070	 InterTube tracheal tube, uncuffed	7.0 mm	10 
9050065	InterTube tracheal tube, uncuffed	6.5 mm	10 
9050060	InterTube tracheal tube, uncuffed	6.0 mm	10 
9050055	InterTube tracheal tube, uncuffed	5.5 mm	10 
9050050	InterTube tracheal tube, uncuffed	5.0 mm	10 
9050045	InterTube tracheal tube, uncuffed	4.5 mm	10 
9050040	InterTube tracheal tube, uncuffed	4.0 mm	10 
9050035	InterTube tracheal tube, uncuffed	3.5 mm	10 
9050030	InterTube tracheal tube, uncuffed	3.0 mm	10 
9050025	InterTube tracheal tube, uncuffed	2.5 mm	10 

 Illustrated  Sterile

For the full range of Airway Management devices and accessories visit
www.intersurgical.com/products/airway-management



InterTube is manufactured by Well Lead Medical Co and distributed by Intersurgical Ltd.

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS6.28 • Issue 2 03.25



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com

**UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • India • Ireland • Italy
 Japan • Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA**