




EcoMask™

The greener, safer anaesthetic face mask



Anaesthesia ▪ Anaesthetic Face Masks



EcoMask™ anaesthetic face masks

The greener, safer anaesthetic face mask

EcoMask is a range of a phthalate-free, non-PVC, single use anatomical anaesthetic face masks which are manufactured from polypropylene and thermoplastic elastomer. EcoMask is a lightweight and clearer version of the traditional, reusable rubber face mask, and is available in a range of seven sizes with distinct colour-coding for easy selection.

The greener, safer choice

EcoMask is manufactured from polypropylene and thermoplastic elastomer, contains no PVC and is phthalate-free.



22F

Non-slip rings


For a better grip and tighter seal

Anatomical shape

For a tailored fit

Soft seal

For a comfortable fit



Transparent shell

For excellent visibility

Distinct colour-coding

For easy selection and integration with other Intersurgical airway products

Code	Description	Box Qty.
7096000	EcoMask, anaesthetic face mask, size 6, extra large adult, red seal, 22F	20
7095000	EcoMask, anaesthetic face mask, size 5, large adult, orange seal, 22F	25
7094000	EcoMask, anaesthetic face mask, size 4, medium adult, green seal, 22F	35
7093000	EcoMask, anaesthetic face mask, size 3, small adult, yellow seal, 22F	35
7092000	EcoMask, anaesthetic face mask, size 2, paediatric, white seal, 22F	25
7091000	EcoMask, anaesthetic face mask, size 1, infant, grey seal, no hook ring, 15M	35
7090000	EcoMask, anaesthetic face mask, size 0, neonate, light blue seal, no hook ring, 15M	40
7089000	Anaesthetic mask hook ring for both QuadraLite™ and EcoMask with 22mm port	150

Make an enquiry



Lower environmental impact product



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
Save energy and paper.
If you must print this information sheet please print it double sided.

IS5.21 • Issue 5 07.19