



nFlow™

Infant nasal CPAP system



Critical Care ▪ Breathing Systems



nFlow™

Designed to care

The smallest patients often need the most care. Safety and comfort are vital to help infants with their earliest steps in life.

A choice for infant CPAP

The Intersurgical nFlow system offers a complete care solution for infant CPAP with established 'gas flip' technology, providing a low 'work of breathing' option as an alternative to ventilator CPAP.


nFlow is compatible with flow drivers designed for use with gas flip technology CPAP, providing slightly higher system pressure associated with this treatment.

The nFlow system

The nFlow system is supplied assembled with generator, heated wire limb and 10 mm adjustable Superset style exhaust tube for efficient gas removal. The complete nFlow breathing system is phthalate free including the small bore tube used for monitoring and gas delivery.




- [Make an enquiry](#)
- [Find out more](#)
- [Watch the video](#)





A choice of patient interfaces and bonnets

Patient interfaces and bonnets are supplied separately to avoid waste. Use the size guide provided to choose the correct option for your patient.




12 bonnets

Available separately in twelve colour-coded sizes. Each bonnet includes 3 fixation pads to secure the tube support and bonnet ties.



Patient interface

Masks and prongs are available separately in four size options. Made from soft, comfortable silicone to ensure optimal patient fit with minimal wastage.



B Bonnets



C Nasal masks



D Nasal prongs





We've made ordering easy. To make a complete system you can choose one product from group A and B, and choose a patient interface from C and/or D. At any time you can order separate items according to your needs.



Code	Description		Box Qty.
A nFlow breathing systems			
470000	nFlow infant CPAP heated wire breathing system. ≥ 1.6 m		12
470001	nFlow infant CPAP heated wire breathing system. ≥ 1.8 m		10
4700016	nFlow infant CPAP heated wire breathing system with detachable monitoring line. ≥ 1.8 m		10
4708000	nFlow silencer kit		50
4700011	nFlow replacement bonnet ties		60
B nFlow bonnets			
		Colour	Size
4907009	nFlow, infant bonnet, navy		Size 9 (40 - 42 cm)
4907008	nFlow, infant bonnet, turquoise		Size 8 (38 - 40 cm)
4907007	nFlow, infant bonnet, orange		Size 7 (36 - 38 cm)
4907006	nFlow, infant bonnet, red		Size 6 (34 - 36 cm)
4907005	nFlow, infant bonnet, green		Size 5 (32 - 34 cm)
4907004	nFlow, infant bonnet, light orange		Size 4 (30 - 32 cm)
4907003	nFlow, infant bonnet, blue		Size 3 (28 - 30 cm)
4907002	nFlow, infant bonnet, yellow		Size 2 (26 - 28 cm)
4907001	nFlow, infant bonnet, dark grey		Size 1 (24 - 26 cm)
4907110	nFlow, infant bonnet, pink		Size 0 (22 - 24 cm)
4907100	nFlow, infant bonnet, grey		Size 00 (20 - 22 cm)
4907000	nFlow, infant bonnet, white		Size 000 (18 - 20 cm)
C nFlow nasal masks			
		Colour	Size
4716000	nFlow nasal mask		Extra large
4706000	nFlow nasal mask		Large
4705000	nFlow nasal mask		Medium
4704000	nFlow nasal mask		Small
D nFlow nasal prongs			
		Colour	Size
4703000	nFlow nasal prong		Large
4702000	nFlow nasal prong		Medium
4701000	nFlow nasal prong		Small
4711000	nFlow nasal prong		Extra small

IS1.16 • issue 9 10.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.